




DESS

> DENTAL SMART SOLUTIONS



Quick reference
guide 2022

PRODUCT LIST



HEALING ABUTMENTS

04



DESSLoc®

17



C-Base®

31



CASTABLES

44



SCAN ABUTMENTS

06



MULTI-UNIT
ABUTMENTS

19



ABUTMENTS

33



TI-BASE CASTABLES

46



TRANSFERS

09



TI-BASES

23




SCREWS

34



AURUM CASTABLES

49



ANALOGUES

12



AURUMBase®

27



PRE-MILLED
BLANKS

39




HOLDERS

51



TEMPORARY
ABUTMENTS

16



ELLIPTIBase®

30



CoCr BASES

42



TOOLS

52

Healing Abutments

Our healing abutments have been manufactured in titanium Grade 5 ELI.
They are identified with the platform and gingival height and have the same connection as the implant screw.
- Manufactured in titanium Grade 5 ELI.



A) Single Crown
B) For Bridges & structures


Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip®	NP (3.0) A	(A) NP (3.5) (B)	(A) RP (4.2) (B)	(A) NP (3.5) (B) WP (5.5)			
		3,0mm	20.068/3	20.039/3	20.047/3	20.040/3	20.048/3	20.069/3	20.070/3
		4,0mm	20.068/4	20.039/4	20.047/4	20.040/4	20.048/4	20.069/4	20.070/4
		5,0mm	20.068/5	20.039/5	20.047/5	20.040/5	20.048/5	20.069/5	20.070/5
	Nobel Replace Select™ Trilobe	UniGrip®	NP (3.5)	RP (4.3)	WP (5.0)				
		3,0mm	20.004/3	20.005/3	20.006/3				
		4,0mm	20.004/4	20.005/4	20.006/4				
	Nobel Brånemark® External Hex Universal	UniGrip®	NP (3.5)	RP (4.1)	WP (5.0)				
		3,0mm	20.001/3	20.002/3	20.003/3				
		4,0mm	20.001/4	20.002/4	20.003/4				
5,0mm		20.001/5	20.002/5	20.003/5					
Multi-Unit® Multiunit Connection	UniGrip®	RP (4.8)	WP (6.5)						
		20.100	20.054						
	Hex 1,20mm	20.100H120							
	Hex 1,22mm	20.100H122							
	Hex 1,27mm	20.100H127							
	Torx® Ball	20.100TB							
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)				
		3,0mm	20.017/3	20.018/3	20.019/3				
		4,0mm	20.017/4	20.018/4	20.019/4				
		5,0mm	20.017/5	20.018/5	20.019/5				
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)				
		3,0mm	20.011/3	20.012/3	20.013/3				
		4,0mm	20.011/4	20.012/4	20.013/4				
	3i Certain® Implants Internal Hex “Click”	5,0mm	20.011/5	20.012/5	20.013/5				
		Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)				
		3,0mm	20.014/3	20.015/3	20.016/3				
	4,0mm	20.014/4	20.015/4	20.016/4					
	5,0mm	20.014/5	20.015/5	20.016/5					
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0			
		3,0mm		20.065/3	20.066/1				
		4,0mm		20.065/4	20.066/4				
		5,0mm		20.065/5	20.066/5				
	Biohorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)			
		3,0mm		20.087/3	20.088/3	20.089/3			
		4,0mm		20.087/4	20.088/4	20.089/4			
	Biohorizons® Implants External External Hex BH	5,0mm		20.087/5	20.088/5	20.089/5			
		Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)				
		3,0mm	20.020/3	20.021/3	20.022/3				
	4,0mm	20.020/4	20.021/4	20.022/4					
	5,0mm	20.020/5	20.021/5	20.022/5					

Dentsply® Sirona®	Astratech Implants System™ EV Conic EVO	Hex ø 1,27mm 3,5mm 4,5mm 6,5mm	EV (3.0) Ø4 20.406/3 Ø4 20.406/4 Ø4 20.406/6	EV (3.6) Ø4 20.059/3 Ø4 20.059/4 Ø4 20.059/6	EV (4.2) Ø5 20.060/3 Ø5 20.060/4 Ø5 20.060/6	EV (4.8) Ø5 20.061/3 Ø6,5 20.061/4 Ø6,5 20.061/6	EV (5.4) Ø5 20.407/3 Ø6,5 20.407/4 Ø6,5 20.407/6			
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm 3,0mm 4,0mm 5,0mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)					
		Anatomic vers. 4,0mm 4,0mm 5,0mm		20.024/3 20.024/4 20.024/5	20.025/3 20.025/4 20.025/5					
				Ø4,5 20.424/4						
				Ø5,5 20.524/4	20.525/4					
					Ø4,5 20.625/4					
		Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm 3,0mm 4,0mm 5,0mm	NP (3.4) 20.036/3 20.036/4 20.036/5	RP (3.8) 20.037/3 20.037/4 20.037/5	WP (4.5) 20.038/3 20.038/4 20.038/5	5.5			
		Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm 1,5mm 2,5mm		20.046/1 20.046/3					
	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm 3,0mm 4,0mm 5,0mm	NP (narrow)	SP (Std) 20.017/3 20.17/4 20.017/5	WP (wide) 20.018/3 20.018/4 20.018/5	Multi-unit 4.8				
	Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx® 2,0mm 3,0mm 4,0mm On SynOcta®	NNC (3.5)	RN (4.8) 20.009/2 20.009/3 20.009/4	WN (6.5) 20.045/2 20.045/3 20.045/4				
					20.010	20.049				
			Straumann® Bone Level Conical BL	Torx® 3,0mm 4,0mm 5,0mm	NP (3.3)	RP (4.1) 20.041/3 20.041/4 20.041/5	NP/RP 20.042/3 20.042/4 20.042/3			
				Straumann® BLX Conical BLX	Torx® 0,8mm ^{A/B} 1,5mm ^A 1,5mm ^B 1,5mm ^{A/B} 2,5mm ^A 2,5mm ^B 3,5mm ^A 3,5mm ^B	RB (3.5)	WB (4.5) 20.463/1 20.461/2 20.462/2 20.463/2 20.461/3 20.462/3 20.461/4 20.462/4	RB/WB		

^A For Single Crown Restoration
^B For Bridges and Structures Restoration

Scan Abutments

Take a digital impression with our desktop, intraoral or Ti-based scanbodies.
Our comprehensive range, combined with our free libraries for ExoCad, 3Shape and Dental Wings, allows you to design with unique precision.



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
			50.068	50.041	50.042	50.069	
			52.068	52.041	52.042	52.069	
			53.100-P10	53.100-P10	53.100-P10	53.100-P10	non-engaging
			53.200-P10	53.200-P10	53.200-P10	53.200-P10	engaging
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
			50.004	50.005	50.006	50.097	
			52.004	52.005	52.006	52.097	
			53.100-P10	53.100-P10	53.100-P10	53.100-P10	non-engaging
			53.200-P10	53.200-P10	53.200-P10	53.200-P10	engaging
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
			50.001	50.002	50.003		
			52.001	52.002	52.003		
			53.100-P10	53.100-P10	53.100-P10		non-engaging
			53.200-P10	53.200-P10	53.200-P10		engaging
	Multi-Unit® Multiunit Connection	UniGrip®	RP (4.8)	WP (6.5)			
		50.007	50.047				
		52.007	52.047				
		53.107	53.154				
		53.107	53.155				
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
			50.017	50.018	50.019		
			52.017	52.018	52.019		
			53.100-P10	53.100-P10	53.100-P10		
			53.200-P10	53.200-P10	53.200-P10		
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
			50.011	50.012	50.013		
			52.011	52.012	52.013		
			53.100-P10	53.100-P10	53.100-P10		
			53.200-P10	53.200-P10	53.200-P10		
	3i Certain® Implants Internal Hex “Click”	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		
			50.014	50.015	50.016		
		52.014	52.015	52.016			
		53.100-P10	53.100-P10	53.100-P10			
		53.200-P10	53.200-P10	53.200-P10			
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.6)	EV (4.2)	EV (4.8)		
			50.059	50.060	50.061		
		On Uniabutment®					
			52.059	52.060	52.061		
			53.100-P10	53.100-P10	53.148-P10		
			53.200-P10	53.200-P10	53.248-P10		

Dentsply® Sirona®	Astra Tech Osseospeed™ Implants Internal Hex Conic	 Hex ø 1,27mm		NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	
				50.023	50.024	50.025	
				52.023	52.024	52.025	
			20°	50.036			
			45°	50.037			
			20°	52.036			
		45°	52.037				
			53.100-P10	53.100-P10	53.148-P10	non-engaging	
			53.200-P10	53.200-P10	53.248-P10	engaging	
	Dentsply Friadent® Xive® Implants Internal Hex FD	 Hex ø 1,22mm		NP (3.4)	RP (3.8)	WP (4.5)	5.5
			50.038	50.039	50.040	50.071	
			52.038	52.039	52.040	52.071	
			53.100-P10	53.100-P10	53.100-P10	53.100-P10	
			53.200-P10	53.200-P10	53.200-P10	53.200-P10	
							
Dentsply Ankylos® C/X Implants Internal ANK	 Hex ø 1,00mm						
						50.049	
						52.049	
						53.100-P10	
					53.200-P10		
MIS® Seven Implants Internal Hex MI	 Hex ø 1,27mm		NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8	
			50.056	50.017	50.018		
			52.056	52.017	52.018	52.085	
			53.100-P10	53.100-P10	53.100-P10	non-engaging	
MIS® C1 Implants Conic MC	 Hex ø 1,27mm		NP (yellow)	SP (purple)	WP (green)	Multi-unit 4.8	
			52.056	52.017	52.018	52.085	
			53.100-P10	53.100-P10	53.100-P10	non-engaging	
			53.200-P10	53.200-P10	53.200-P10	engaging	
MIS® V3 Implants Conic MV	 Hex ø 1,27mm		NP (blue)	SP (purple)	Multi-unit 4.8		
			52.084	52.082	52.085		
			53.100-P10	53.100-P10		non-engaging	
		53.200-P10	53.200-P10		engaging		
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	 Torx®		NNC (3.5)	RN (4.8)	WN (6.5)	
				50.070	50.009	50.045	
		 On Abutment			50.010	50.048	
				52.070	52.009	52.045	
		 On Abutment			52.010	52.048	
				53.100-P10	53.110-P10	53.148-P10	non-engaging
	Straumann® Bone Level Conical BL			53.200-P10	53.210-P10	53.248-P10	engaging
		 Torx®		NP (3.3)	RP (4.1)	NP/RP	
				50.043	50.044		
		 On Abutment SRA				50.098	
				52.043	52.044		
		 On Abutment SRA				52.098	
	Straumann® BLX Conical BLX			53.100-P10	53.100-P10		non-engaging
				53.200-P10	53.200-P10		engaging
		 Torx®		RB (3.5)	WB (4.5)		
				52.062	52.062		
Neodent Grand Morse® NEO GM			53.100-P10	53.100-P10		non-engaging	
			53.200-P10	53.200-P10		engaging	
					GM		
				52.062			
				53.100-P10	non-engaging		
				53.200-P10	engaging		

Scan Abutments

Straumann® Group	Anthogyr Axiom® Implants Conic Axio		52.062				non-engaging				
			53.100-P10								
			53.200-P10					engaging			
Henry Schein® Group	Camlog® Implants Internal Cam		3.3				3.8	4.3	5.0	non-engaging	
			50.064				50.065	50.066	50.067		
			52.064				52.065	52.066	52.067		
			53.100-P10				53.100-P10	53.100-P10	53.100-P10		engaging
	Conelog® Implants Internal Con		NP (3.3)				RP (3.8)	WP (5.0)	engaging		
			52.074				52.075	52.076			
			53.100-P10				53.100-P10	53.100-P10			
			53.200-P10				53.200-P10	53.200-P10			
	BioHorizons® Implants Internal Internal Hex BH		3.0				NP (3.5)	RP (4.5)	WP (5.7)	non-engaging	
			52.099				52.087	52.088	52.089		
							53.100-P10	53.100-P10	53.100-P10		engaging
							53.200-P10	53.200-P10	53.200-P10		
	BioHorizons® Implants External External Hex BH		NP (3.7)				RP (4.2)	WP (5.2)	non-engaging		
			50.020				50.021	50.022			
			52.020				52.021	52.022			
			53.100-P10				53.100-P10	53.100-P10		engaging	
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss		MINI				STD	non-engaging			
			50.072				50.073				
			52.072				52.073				
Megagen®	Megagen® Anyridge Implants Conic Anyr		53.100-P10				53.100-P10	engaging			
			53.200-P10				53.200-P10				
Medentis®	Medentis ICX® Implants Conic IC		50.080				50.080	non-engaging			
			52.080				52.080				
			53.100-P10				53.100-P10				
GlobalD®	GlobalD® Implants Conic Ink		53.200-P10				53.200-P10	engaging			
Dentium®	Superline™ & Implantium® Conic Dent		52.062				52.062	non-engaging			
			53.100-P10				53.100-P10				
			53.200-P10				53.200-P10		engaging		

Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27	3.25/3.75	4.1	4.5	5.5	non-engaging engaging	
			52.091	52.092	52.093	52.094		
			53.100-P10	53.100-P10	53.100-P10	53.100-P10		
			53.200-P10	53.200-P10	53.200-P10	53.200-P10		
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22	NP		RP		non-engaging engaging	
			52.078					
			53.100-P10					
			53.200-P10					
BTI®	BTI® Implants Internal BT	Hex ø 1,22	4.1				non-engaging engaging	
			52.079					
			53.100-P10					
			53.200-P10					
Phibo®	Phibo® TSH® External Hex PH		NP	SP	WP	non-engaging engaging		
			52.030	52.031	52.031			
			53.100-P10	53.100-P10	53.100-P10			
			53.200-P10	53.200-P10	53.200-P10			

Transfers

With a body manufactured in stainless steel AISI 303, and the open and closed tray screws included manufactured in titanium, the DESS Transfer has been designed with a flat face to improve retention in the model and grooves to easily splint them together. The Transfer screw can be tightened using different method-specific system screwdriver, spatula, sharpened object or manually.



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)
			17.068	17.044	17.045	17.069
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0
			17.004	17.005	17.006	17.097
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)	
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA					
			17.017	17.018	17.019	
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)	
			17.011	17.012	17.013	

Transfers

Zimmer® 3i® Biomet®	3i Certain® Implants Internal Hex “Click”	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)	
			17.014	17.015	17.016	
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.6)	EV (4.2)	EV (4.8)
			17.058	17.059	17.060	17.061
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm		NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)
		20°	17.036	17.023	17.024	17.025
		45°	17.037			
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5
			17.038	17.039	17.040	17.071
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm				
		GH 1,5mm				17.049
		GH 2,5mm				17.250
	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	
			17.056	17.017	17.018	
	MIS® C1 Implants Conic MC	Hex ø 1,27mm	NP (yellow)	SP (purple)	WP (green)	Multi-unit 4.8
			17.081	17.082	17.083	
						17.085
	MIS® V3 Implants Conic MV	Hex ø 1,27mm	NP (blue)	SP (purple)		
		17.084	17.082			
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
			17.070	17.009	17.053	
		On Abutment		17.010	17.048	
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)		
			17.046	17.047		
	Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)		
			17.046	17.046		

Straumann® Group	Neodent Grand Morse® NEO GM	NEO				GM
						17.035
	Anthogyr Axiom® Implants Conic Axio	Hex ø 1,18mm				
						17.077
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0
			17.064	17.065	17.066	17.067
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)	
			17.074	17.075	17.076	
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
				17.087	17.088	17.089
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD		
			17.072	17.073		
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm				
						17.057
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm				
						17.057
GlobalID®	GlobalID® Implants Conic Ink	Hex ø 1,22mm				
						17.086
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28				
						17.090


Transfers

Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27	3.25/3.75	4.1	4.5	5.5
			17.091	17.092	17.093	17.094
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22		NP	RP	
			17.078			
BTI®	BTI® Implants Internal BT	Hex ø 1,22	4.1			
			17.079			
Phibo®	Phibo® TSH® External Hex PH		NP	SP	WP	
			17.030	17.031	17.032	
DESSLoc®	DESSLoc®					17.054

Analogue

Manufactured in surgical stainless steel (metal finish) or titanium Grade 5 ELI colour coded (anodised) to match the implant platform.

- Self-locking
- Excellent stability in stone models
- Digital analogue includes a screw to secure the analogue to the bottom of the model



Traditional: xx.xxx
Digital: xx.xxx/D

Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex		3.0	NP (3.5)	RP (4.3)	WP (5.5)
			14.068/D	14.041/D	14.042/D	14.069/D
			14.068	14.041	14.042	14.069
				14.041-P10	14.042-P10	14.069-P10
	Nobel Replace Select™ Trilobe		NP (3.5)	RP (4.3)	WP (5.0)	6.0
			14.004/D	14.005/D	14.006/D	14.097/D
			14.004	14.005	14.006	14.097
			14.004-P10	14.005-P10		
	Nobel Brånemark® External Hex Universal		NP (3.5)	RP (4.1)	WP (5.0)	
			14.001/D	14.002/D	14.003/D	
			14.001	14.002	14.003	
				14.002-P10		
	Multi-Unit® Multiunit Connection		RP (4.8)	WP (6.5)		
			14.007/D	14.047/D		
			14.007	14.047		
			14.007-P10	14.047-P10		

Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex		3.0	NP (3.5)	RP (4.3)	WP (5.5)
			14.068/D	14.041/D	14.042/D	14.069/D
			14.068	14.041	14.042	14.069
				14.041-P10	14.042-P10	14.069-P10
	Nobel Replace Select™ Trilobe		NP (3.5)	RP (4.3)	WP (5.0)	6.0
			14.004/D	14.005/D	14.006/D	14.097/D
			14.004	14.005	14.006	14.097
			14.004-P10	14.005-P10		
	Nobel Brånemark® External Hex Universal		NP (3.5)	RP (4.1)	WP (5.0)	
			14.001/D	14.002/D	14.003/D	
			14.001	14.002	14.003	
				14.002-P10		
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA		NP (3.5)	RP (4.5)	WP (5.7)	
			14.017/D	14.018/D	14.019/D	
			14.017	14.018	14.019	
			14.017-P10	14.018-P10		
	3i Osseotite® External Hex USA		NP (3.4)	RP (4.1)	WP (5.0)	
			14.011/D	14.012/D	14.013/D	
			14.011	14.012	14.013	
				14.012-P10		
	3i Certain® Implants Internal Hex “Click”		NP (3.4)	RP (4.1)	WP (5.0)	
			14.014/D	14.015/D	14.016/D	
			14.014	14.015	14.016	
				14.015-P10		
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Digital 1un.	Of Uniabutment®	EV (3.6)	EV (4.2)	EV (4.8)
				14.059/D	14.060/D	14.061/D
		Standard 1un.		14.058/D		
				14.058	14.059	14.060
						14.061
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Digital 1un.	Of Uniabutment®	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)
		Digital 1un. 20°		14.023/D	14.024/D	14.025/D
		Digital 1un. 45°		14.037/D		
		Std. 1un.			14.023	14.024
		Std. 10un.			14.024	14.025
		Std. 1un. 20°		14.023-P10	14.024-P10	14.025-P10
		Std. 10un. 20°		14.036		
		Std. 10un. 20°		14.036-P10		
		Std. 1un. 45°		14.037		
		Std. 10un. 45°		14.037-P10		
	Dentsply Friadent® Xive® Implants Internal Hex FD	Digital 1un.		NP (3.4)	RP (3.8)	WP (4.5)
		Std. 1un.		14.038/D	14.039/D	14.071/D
		Std. 10un.		14.038	14.039	14.040
				14.038-P10	14.039-P10	
	Dentsply Ankylos® C/X Implants Internal ANK	Digital 1un. ø 3,5				14.048/D
		Digital 1un. ø 4,5				14.049/D
		Std. 1un. ø 3,5				14.048
		Std. 10un. ø 3,5				14.048-P10
		Std. 1un. ø 4,5				14.049
		Std. 10un. ø 4,5				14.049-P10

Analogues


Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	Digital 1un.	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
		Digital 1un.	14.056/D	14.095/D	14.096/D	
		Std. 1un.	14.056	14.095	14.096	14.085/D mis Multi-Unit
		Std. 10un.	14.056-P10	14.095-P10	14.096-P10	
	MIS® C1 Implants Conic MC		NP (yellow)	SP (purple)	WP (green)	Multi-unit 4.8
			14.081/D	14.082/D	14.083/D	
			14.081	14.082	14.083	14.085/D
						14.085
	MIS® V3 Implants Conic MV		NP (blue)	SP (purple)	Multi-unit 4.8	
			14.084/D	14.082/D		
					14.085/D	
			14.084	14.082	14.085	
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
			14.070/D	14.009/D	14.045/D	
		On Abutment		14.010/D	14.046/D	
			14.070	14.009	14.045	
	Straumann® Bone Level Conical BL			14.009-P10		
				14.010	14.046	
				14.010-P10		
	Straumann® BLX Conical BLX	Torx®	NP (3.3)	RP (4.1)	NP/RP	
		On Abutment SRA	14.043/D	14.044/D		
					14.098/D	
		On Abutment SRA	14.043	14.044		
	Straumann® BLX Conical BLX		14.043-P10	14.044-P10		
					14.098	
	Neodent Grand Morse® NEO GM		RB (3.5)	WB (4.5)		
			14.062/D	14.063/D		
			14.062	14.063		
	Neodent Grand Morse® NEO GM	NEO			GM	
					14.035/D	
					14.035	
	Anthogyr Axiom® Implants Conic Axio				14.077/D	
					14.077	
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0
			14.064/D	14.065/D	14.066/D	14.067/D
			14.064	14.065	14.066	14.067
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)	
			14.074/D	14.075/D	14.076/D	
			14.074	14.075	14.076	

Henry Schein® Group	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
			14.099/D	14.087/D	14.088/D	14.089/D
	BioHorizons® Implants External External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)	
			14.020/D	14.021/D	14.022/D	
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm				
			14.072/D	14.073/D		
			14.072	14.073		
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm				
			14.057/D			
			14.057			
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm				
			14.080/D			
			14.080			
GlobalD®	GlobalD® Implants Conic Ink	Hex ø 1,22mm				
			14.086/D			
			14.086			
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28				
			14.090/D			
			14.090			
Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27	3.25/3.75	4.1	4.5	5.5
			14.091/D	14.092/D	14.093/D	14.094/D
			14.091	14.092	14.093	14.094
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22		NP	RP	
			14.078/D			
			14.078			
BTI®	BTI® Implants Internal BT	Hex ø 1,22	4.1			
			14.079/D			
			14.079			
Phibo®	Phibo® TSH® External Hex PH		NP	SP	WP	
			14.030/D	14.031/D	14.031/D	
			14.030	14.031	14.031	

Temporary Abutments

Manufactured in titanium Grade 5 ELI. Our Temporary abutment comes with a 10,5 mm chimney height that can be cut to the required length. With parallel stretch marks to improve adhesion of the material.

- Screw included




Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
				24.041	24.042		non-engaging
				25.041	25.042		engaging
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
			24.004	24.005	24.006		non-engaging
			25.004	25.005	25.006		engaging
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
			24.001	24.002	24.003		non-engaging
			25.001	25.002	25.003		engaging
	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			
			24.007	24.030			non-engaging
			24.007H120				non-engaging
			24.007H122				non-engaging
			24.007H127				non-engaging
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
			24.017	24.018	24.019		non-engaging
			25.017	25.018	25.019		engaging
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
			24.011	24.012	24.013		non-engaging
			25.011	25.012	25.013		engaging
	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		
			24.014	24.015	24.016		non-engaging
			25.014	25.015	25.016		engaging
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.6)	EV (4.2)	EV (4.8)	
				24.059	24.060	24.061	non-engaging
				25.059	25.060	25.061	engaging
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm					
			NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)		
				24.024	24.025		non-engaging
				25.024	25.025		engaging
				24.028			
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm					
			NP (3.4)	RP (3.8)	WP (4.5)	5.5	
			24.038	24.039	24.040		non-engaging

Dentsply® Sirona®	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm				24.048	non-engaging
						24.250	non-engaging
						25.048	engaging
						25.250	engaging
	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8	
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®					
	Straumann® Bone Level Conical BL	On Abutment On Abutment Torx®					
		On Abutment SRA On Abutment SRA					
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm					
	BioHorizons® Implants External External Hex BH	Hex ø 1,27mm					
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm					

DESSLoc®

Removable overdenture attachment. With a zirconium nitride coating that produces a stronger abutment and reduces inflammation and plaque adhesion. Available in up to 5 different heights.



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
		GH 1,0mm		88.039/1	88.040/1		
		GH 2,0mm		88.039/2	88.040/2		
		GH 3,0mm		88.039/3	88.040/3		
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		GH 1,0mm	88.004/1	88.005/1			
		GH 2,0mm	88.004/2	88.005/2			
		GH 3,0mm	88.004/3	88.005/3			
		GH 4,0mm	88.004/4	88.005/4			

DESSLoc®

Nobel Biocare®	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)
		GH 1,0mm	88.001/1	88.002/1	
		GH 2,0mm	88.001/2	88.002/2	
		GH 3,0mm	88.001/3	88.002/3	
		GH 4,0mm		88.002/4	
		GH 5,0mm		88.002/5	
	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8) 88.007	WP (6.5)	
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)
		GH 1,0mm	88.017/1	88.018/1	88.019/1
		GH 2,5mm	88.017/2	88.018/2	88.019/2
		GH 3,5mm	88.017/3	88.018/3	88.019/3
		GH 4,5mm	88.017/4	88.018/4	88.019/4
		GH 5,5mm	88.017/5	88.018/5	
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)
		GH 1,0mm	88.011/1	88.002/1	88.013/1
		GH 2,0mm	88.011/2	88.002/2	88.013/2
		GH 3,0mm	88.011/3	88.002/3	88.013/3
		GH 4,0mm	88.011/4	88.002/4	88.013/4
		GH 5,0mm	88.011/5	88.002/5	
	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)
		GH 1,0mm	88.014/1	88.015/1	88.016/1
		GH 2,0mm	88.014/2	88.015/2	88.016/2
		GH 3,0mm	88.014/3	88.015/3	88.016/3
		GH 4,0mm	88.014/4	88.015/4	88.016/4
		GH 5,0mm	88.014/5	88.015/5	
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.6)	EV (4.2)	EV (4.8)
		GH 1,0mm	88.059/1	88.060/1	88.061/1
		GH 2,0mm	88.059/2	88.060/2	88.061/2
		GH 3,0mm	88.059/3	88.060/3	88.061/3
		GH 4,0mm	88.059/4	88.060/4	88.061/4
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)
		GH 1,0mm	88.023/1	88.024/1	88.025/1
		GH 2,0mm	88.023/2	88.024/2	88.025/2
		GH 3,0mm	88.023/3	88.024/3	88.025/3
		GH 4,0mm	88.023/4	88.024/4	88.025/4
		GH 5,0mm			88.025/5
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)
		GH 1,0mm	88.036/1	88.037/1	5.5
		GH 2,0mm	88.036/2	88.037/2	
		GH 3,0mm	88.036/3	88.037/3	
		GH 4,0mm	88.036/4	88.037/4	
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm			
		GH 2,0mm		88.046/2	
		GH 3,0mm		88.046/3	
		GH 4,0mm		88.046/4	

Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
		GH 1,0mm	88.017/1	88.018/1		
		GH 2,5mm	88.017/2	88.018/2		
		GH 3,5mm	88.017/3	88.018/3		
		GH 4,5mm	88.017/4	88.018/4		
		GH 5,5mm	88.017/5			
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
		GH 1,0mm		88.009/1		
		GH 2,0mm		88.009/2		
		GH 3,0mm		88.009/3		
		GH 4,0mm		88.009/4		
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP	
		GH 1,0mm		88.042/1		
		GH 2,0mm	88.041/2	88.042/2		
		GH 3,0mm	88.041/3	88.042/3		
		GH 4,0mm	88.041/4	88.042/4		
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0
		GH 1,0mm		88.065/1	88.066/1	
		GH 2,0mm		88.065/2	88.066/2	
		GH 3,0mm		88.065/3	88.066/3	
		GH 4,0mm		88.065/4	88.066/4	
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
		GH 1,0mm		88.087/1	88.088/1	
		GH 2,5mm		88.087/2	88.088/2	
		GH 3,5mm		88.087/3	88.088/3	
		GH 4,5mm		88.087/4	88.088/4	
		GH 5,5mm		88.087/5		

Multi-Units

Available in 14 compatibilities, different gingival heights, and a choice of 0°, 17° and 30° angulation, the DESS® range offers you the most complete range in the market.



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
		GH 1,2mm 0°		40.007/1	40.008/1		non-engaging
		GH 2,5mm 0°		40.007/2	40.008/2		non-engaging
		GH 3,5mm 0°		40.007/3	40.008/3		non-engaging
		GH 4,5mm 0°			40.008/4		non-engaging
		GH 2,5mm 17°		41.707/2	41.708/2		engaging
		GH 3,5mm 17°		41.707/3	41.708/3		engaging
		GH 3,5mm 30°		43.007/3	43.008/3		engaging
		GH 4,5mm 30°		43.007/4	43.008/4		engaging

Multi-Units

Nobel Biocare®	Nobel Replace Select™ Trilobe	UniGrip		NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		GH 1,2mm	0°	40.004/1	40.005/1	40.006/1		non-engaging
		GH 2,5mm	0°	40.004/2	40.005/2	40.006/2		non-engaging
		GH 3,5mm	0°	40.004/3	40.005/3	40.006/3		non-engaging
		GH 4,5mm	0°		40.005/4			non-engaging
		GH 4,5mm	0°		40.005/5			non-engaging
		GH 2,5mm	17°	41.704/2	41.705/2			engaging
		GH 3,5mm	17°	41.704/3	41.705/3			engaging
		GH 2,5mm	17°		41.705/4			engaging
		GH 3,5mm	30°		43.005/4			engaging
		GH 3,5mm	30°		43.005/5			engaging
		GH 4,5mm	30°		43.005/5R			engaging
	Nobel Brånemark® External Hex Universal	UniGrip		NP (3.5)	RP (4.1)	WP (5.0)		
		GH 1,0mm	0°	40.001/1	40.002/1	40.003/1		non-engaging
		GH 2,0mm	0°	40.001/2	40.002/2	40.003/2		non-engaging
		GH 3,0mm	0°	40.001/3	40.002/4	40.003/3		non-engaging
		GH 4,0mm	0°		40.002/4			non-engaging
		GH 5,0mm	0°		41.002/5			engaging
		GH 2,0mm	17°	41.701/2	41.702/2			engaging
		GH 3,0mm	17°	41.701/3	41.702/3			engaging
		GH 4,0mm	17°		43.702/4			engaging
		GH 4,0mm	30°		43.002/4			engaging
		GH 5,0mm	30°		43.002/5			engaging
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm		NP (3.5)	RP (4.5)	WP (5.7)		
		GH 1,0mm	0°		40.019/1	40.020/1		non-engaging
		GH 2,0mm	0°		40.019/2	40.020/2		non-engaging
		GH 3,0mm	0°		40.019/3	40.020/3		non-engaging
		GH 4,0mm	0°		40.019/4			non-engaging
		GH 5,0mm	0°		40.019/5			non-engaging
		GH 2,5mm	17°		41.719/2	41.720/2		engaging
		GH 3,5mm	17°		41.719/3	41.720/3		engaging
		GH 2,5mm	30°		43.019/3	43.020/3		engaging
		GH 3,5mm	30°		43.019/4	43.020/4		engaging
	3i Osseotite® External Hex USA	Hex ø 1,27mm		NP (3.4)	RP (4.1)	WP (5.0)		
		GH 1,0mm	0°	40.009/1	40.017/1			non-engaging
		GH 2,0mm	0°	40.009/2	40.017/2			non-engaging
		GH 3,0mm	0°	40.009/3	40.017/3			non-engaging
		GH 4,0mm	0°		40.017/4			non-engaging
		GH 5,0mm	0°		40.017/5			non-engaging
		GH 2,0mm	17°		41.712/2			engaging
		GH 3,0mm	17°		41.712/3			engaging
		GH 4,0mm	17°		41.712/4			engaging
		GH 2,0mm	30°		43.012/2			engaging
		GH 3,0mm	30°		43.012/3			engaging
		GH 4,0mm	30°		43.012/4			engaging

Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.6)	EV (4.2)	EV (4.8)	
		GH 1,0mm		40.059/1	40.060/1	40.061/1	non-engaging
		GH 2,0mm		40.059/2	40.060/2	40.061/2	non-engaging
		GH 3,0mm		40.059/3	40.060/3	40.061/3	non-engaging
		GH 5,0mm		40.059/5	40.060/5	40.061/5	non-engaging
		Hex ø 1,27mm		EV (3.6)	EV (4.2)	EV (4.8)	
		GH 2,5mm 0°		40.027/1	40.036/1	40.037/1	non-engaging
		GH 3,5mm 0°		40.027/2	40.036/2	40.037/2	non-engaging
		GH 4,5mm 0°		40.027/3	40.036/3	40.037/3	non-engaging
		GH 2,5mm 17°		41.727/3	41.736/3	41.737/3	engaging
		GH 3,5mm 17°		41.727/4	41.736/4	41.737/4	engaging
		GH 3,5mm 30°		43.027/4	43.036/4	43.037/4	engaging
		GH 4,5mm 30°		43.027/5	43.036/5	43.037/5	engaging
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm		NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	
		GH 0,5mm 20°			40.013/0	40.015/0	non-engaging
		GH 1,0mm 20°			40.013/1	40.015/1	non-engaging
		GH 2,0mm 20°			40.013/2	40.015/2	non-engaging
		GH 4,0mm 20°			40.013/4	40.015/4	non-engaging
		GH 6,0mm 20°			40.013/6	40.015/6	non-engaging
		GH 0,5mm 45°			40.014/0	40.016/0	non-engaging
		GH 1,0mm 45°			40.014/1	40.016/1	non-engaging
		GH 2,0mm 45°			40.014/2	40.016/2	non-engaging
		GH 4,0mm 45°			40.014/4	40.016/4	non-engaging
		GH 6,0mm 45°			40.014/6	40.016/6	non-engaging
		Hex ø 1,27mm		NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	
		GH 2,5mm 0°			40.024/1	40.025/1	non-engaging
		GH 3,5mm 0°			40.024/2	40.025/2	non-engaging
		GH 4,5mm 0°			40.024/3	40.025/3	non-engaging
		GH 2,5mm 17°			41.724/3	41.725/3	engaging
		GH 3,5mm 17°			41.724/4	41.725/4	engaging
		GH 3,5mm 30°			43.024/4	43.025/4	engaging
		GH 4,5mm 30°			43.024/5	43.025/5	engaging
	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm		NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
		GH 1,0mm 0°			40.019/1	40.020/1	non-engaging
		GH 2,0mm 0°			40.019/2	40.020/2	non-engaging
		GH 3,0mm 0°			40.019/3	40.020/3	non-engaging
		GH 4,0mm 0°			40.019/4		non-engaging
		GH 5,0mm 0°			40.019/5		non-engaging
		GH 2,5mm 17°			41.719/2	41.720/2	engaging
		GH 3,5mm 17°			41.719/3	41.720/3	engaging
		GH 3,5mm 30°			43.019/3	43.020/3	engaging
		GH 4,5mm 30°			43.019/4	43.020/4	engaging
Straumann® Group	Straumann® Bone Level Conical BL	Torx®		NP (3.3)	RP (4.1)	NP/RP	
		GH 1,5mm 0°			40.043/1	40.044/1	non-engaging
		GH 2,5mm 0°			40.043/2	40.044/2	non-engaging
		GH 3,5mm 0°			40.043/3	40.044/3	non-engaging
		GH 4,5mm 0°				40.044/4	non-engaging
		GH 2,5mm 17°			41.743/2	41.744/2	engaging
		GH 3,5mm 17°			41.7043/3	41.744/3	engaging
		GH 3,5mm 30°			43.043/3	43.044/3	engaging
		GH 4,5mm 30°			43.043/4	43.044/4	engaging

Multi-Units

Straumann® Group	Straumann® BLX Conical BLX	Torx®		RB (3.5)	WB (4.5)	RB/WB	
		GH 1,5mm	0°		40.062/1		non-engaging
		GH 2,5mm	0°		40.062/2		non-engaging
		GH 3,5mm	0°		40.062/3		non-engaging
		GH 2,5mm	17°		41.762/2		engaging
		GH 3,5mm	17°		41.762/3		engaging
		GH 3,5mm	30°		43.062/3		engaging
		GH 4,5mm	30°		43.062/4		engaging
	Neodent Grand Morse® NEO GM	NEO		GM			
		GH 1,5mm	0°	40.035/1			non-engaging
		GH 2,5mm	0°	40.035/2			non-engaging
		GH 3,5mm	0°	40.035/3			non-engaging
		GH 2,5mm	17°	41.735/2			engaging
		GH 3,5mm	17°	41.735/3			engaging
		GH 3,5mm	30°	43.035/3			engaging
		GH 4,5mm	30°	43.035/4			engaging
Henry Schein® Group	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm		3.0	NP (3.5)	RP (4.5)	WP (5.7)
		GH 1,0mm	0°		40.087/1	40.088/1	non-engaging
		GH 2,0mm	0°		40.087/2	40.088/2	non-engaging
		GH 3,0mm	0°		40.087/3	40.088/3	non-engaging
		GH 4,0mm	0°		40.087/4		non-engaging
		GH 5,0mm	0°		40.087/5		non-engaging
		GH 2,5mm	17°		41.787/2	41.720/2	engaging
		GH 3,5mm	17°		41.787/3	41.720/3	engaging
		GH 3,5mm	30°		43.087/3	43.020/3	engaging
		GH 4,5mm	30°		43.087/4	43.020/4	engaging
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm		MINI	STD		
		GH 1,5mm	0°	40.072/1	40.073/1		non-engaging
		GH 2,5mm	0°	40.072/2	40.073/2		non-engaging
		GH 3,5mm	0°	40.072/3	40.073/3		non-engaging
		GH 2,5mm	17°	41.772/2	41.773/2		engaging
		GH 3,5mm	17°	41.772/3	41.773/3		engaging
		GH 3,5mm	30°	43.072/3	43.073/3		engaging
		GH 4,5mm	30°	43.072/4	43.073/4		engaging
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28					
		GH 1,5mm	0°	40.090/1			non-engaging
		GH 2,5mm	0°	40.090/2			non-engaging
		GH 3,5mm	0°	40.090/3			non-engaging
		GH 2,5mm	17°	41.790/2			engaging
		GH 3,5mm	17°	41.790/3			engaging
		GH 2,5mm	30°	43.090/2			engaging
		GH 3,5mm	30°	43.090/3			engaging

Ti-Bases

Zimmer ® 3i ® Biomet®	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
		GH 0,3mm	15.011	15.012	15.013		non-engaging
		GH 1,5mm	15.211	15.212	15.213		non-engaging
		GH 3,0mm	15.311	15.312	15.313		non-engaging
		GH 0,3mm	16.011	16.012	16.013		engaging
		GH 1,5mm	16.211	16.212	16.213		engaging
		GH 3,0mm	16.311	16.312	16.313		engaging
	3i Certain® Implants Internal Hex “Click”	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		
		GH 0,3mm	15.014	15.015	15.016		non-engaging
		GH 1,5mm	15.214	15.215	15.216		non-engaging
		GH 3,0mm	15.314	15.315	15.316		non-engaging
		GH 0,3mm	16.014	16.015	16.016		engaging
		GH 1,5mm	16.211	16.215	16.216		engaging
		GH 3,0mm	16.314	16.315	16.316		engaging

Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.6)	EV (4.2)	EV (4.8)	EV (5.4)	
		GH 1,0mm		15.059	15.060	15.061	15.407	non-engaging
		GH 3,0mm		15.359	15.360	15.361		non-engaging
		GH 1,0mm		16.059	16.060	16.061	16.407	engaging
		GH 3,0mm		16.359	16.360	16.361		engaging
		On Uniabutment®	15.058					engaging
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)			
		GH 0,9mm			15.025			non-engaging
		GH 1,0mm	15.023					non-engaging
		GH 1,1mm		15.024				non-engaging
		GH 2,0mm	15.023/2	15.024/2	15.025/2			non-engaging
		GH 3,0mm	15.023/3	15.024/3	15.025/3			non-engaging
		GH 0,9mm			16.025			engaging
		GH 1,0mm	16.023					engaging
		GH 1,1mm		16.024				engaging
		GH 2,0mm	16.023/2	16.024/2	16.025/2			engaging
		GH 3,0mm	16.023/3	16.024/3	16.025/3			engaging
		Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5	
	GH 0,3mm		15.038	15.039	15.040	15.071		non-engaging
	GH 1,5mm		15.238	15.239	15.240	15.271		non-engaging
	GH 3,0mm		15.338	15.339	15.340	15.371		non-engaging
	GH 0,3mm		16.038	16.039	16.040	16.071		engaging
	GH 1,5mm		16.238	16.239	16.240	16.271		engaging
	GH 3,0mm		16.338	16.339	16.340	16.371		engaging
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm						
		GH 1,5mm				15.050		non-engaging
		GH 2,5mm				15.250		non-engaging
		GH 1,5mm				16.050		engaging
		GH 2,5mm				16.250		engaging

Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8	
		GH 0,3mm			15.096	15.085	non-engaging
		GH 0,8mm	15.056	15.095			non-engaging
		GH 1,5mm	15.256	15.295	15.296		non-engaging
		GH 3,0mm	15.356	15.395	15.096		non-engaging
		GH 0,3mm			16.096		engaging
		GH 0,8mm	16.056	16.095			engaging
	MIS® C1 Implants Conic MC	Hex ø 1,27mm	NP (yellow)	SP (purple)	WP (green)	Multi-unit 4.8	
		GH 0,3mm				15.085	non-engaging
		GH 1,0mm	15.081	15.082	15.083		non-engaging
		GH 0,3mm				16.085	engaging
	MIS® V3 Implants Conic MV	Hex ø 1,27mm		NP (blue)	SP (purple)	Multi-unit 4.8	
		GH 0,3mm				15.085	non-engaging
		GH 1,0mm		15.084	15.082		non-engaging
		GH 0,3mm					
		GH 1,0mm		16.084	16.082		engaging
		GH 1,0mm					
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)		
		GH 0,3mm			15.045		non-engaging
		GH 0,4mm		15.046			non-engaging
		GH 0,5mm	15.070				non-engaging
		GH 0,3mm			16.045		engaging
	Straumann® Bone Level Conical BL	GH 0,4mm		16.046			engaging
		GH 0,5mm	16.070				
		On Abutment		15.010	15.048		non-engaging
		On Abutment		16.010	16.048		engaging
		Torx®	NP (3.3)	RP (4.1)	NP/RP		
		GH 1,2mm	15.043	15.044			non-engaging
		GH 3,0mm	15.343	15.344			non-engaging
		GH 1,2mm	16.043	16.044			engaging
	Straumann® BLX Conical BLX	GH 3,0mm	16.343	16.344			engaging
		On SRA			15.098		non-engaging
		On SRA			16.098		engaging
		Torx®	RB (3.5)	WB (4.5)	RB/WB		
		GH 0,8mm		15.063			non-engaging
		GH 1,5mm	15.062/2	15.063/3	15.062/2		non-engaging
		GH 2,5mm	15.062/3		15.062/3		non-engaging
		GH 3,5mm	15.062/4		15.062/4		non-engaging
	Neodent Grand Morse® NEO GM	GH 0,8mm		16.063			engaging
		GH 1,5mm	16.062/2	16.063/3	16.062/2		engaging
		GH 2,5mm	16.062/3		16.062/3		engaging
		GH 3,5mm	16.062/4		16.062/4		engaging
		NEO	GM				
		GH 0,8mm	15.035				non-engaging
		GH 1,5mm	15.235				non-engaging
	Anthogyr Axiom® Implants Conic Axio	GH 2,5mm	15.335				non-engaging
		GH 0,8mm	16.035				engaging
		GH 1,5mm	16.235				engaging
		GH 2,5mm	16.335				engaging
		Hex ø 1,18mm					
		GH 2,9mm			15.077		non-engaging
		GH 2,9mm			16.677		engaging


Ti-Bases

Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0	
		GH 0,3mm		15.065	15.066	15.067	non-engaging
		GH 0,5mm	15.064				non-engaging
		GH 1,5mm		15.265	15.266		non-engaging
		GH 3,0mm		15.065	15.066		non-engaging
		GH 0,3mm		16.065	16.066	16.067	engaging
		GH 0,5mm	16.064				engaging
		GH 1,5mm		16.265	16.265		engaging
		GH 3,0mm		16.365	16.366		engaging
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)		
		GH 1,4mm	15.074	15.075	15.076		non-engaging
		GH 3,0mm	15.074/3	15.075/3	15.076/3		non-engaging
		GH 1,4mm	16.074	16.075	16.076		engaging
		GH 3,0mm	16.074/3	16.075/3	16.076/3		engaging
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
		GH 0,3mm			15.088	15.089	non-engaging
		GH 0,8mm		15.087			non-engaging
		GH 0,3mm			16.088	16.089	engaging
	BioHorizons® Implants External External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)		
		GH 0,3mm	15.051	15.052	15.053		non-engaging
		GH 0,3mm	16.051	16.052	16.053		engaging
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD			
		GH 1,0mm	15.072	15.073			non-engaging
		GH 3,0mm	15.372	15.373			non-engaging
		GH 1,0mm	16.072	16.073			engaging
		GH 3,0mm	16.372	16.373			
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm					
		GH 2,0mm			15.057		non-engaging
		GH 3,0mm			15.357		non-engaging
		GH 2,0mm			16.057		engaging
		GH 3,0mm			16.357		engaging
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm			15.080		non-engaging
					16.080		engaging
GlobalID®	GlobalID® Implants Conic Ink				15.086		non-engaging
					16.086		engaging
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28					
		GH 0,9mm		15.090			non-engaging
		GH 1,5mm		15.190			non-engaging
		GH 3,0mm		15.390			non-engaging
		GH 0,9mm		16.090			engaging
		GH 1,5mm		16.190			engaging
		GH 3,0mm		16.390			engaging

Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27	3.25/3.75	4.1	4.5	5.5	
		GH 0,7mm			15.091	15.094	non-engaging
		GH 0,9mm		15.091			non-engaging
		GH 1,1mm	15.091				non-engaging
		GH 0,7mm			16.091	16.094	engaging
		GH 0,9mm		16.091			engaging
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22		NP	RP		
			15.078				non-engaging
			16.078				engaging
BTI®	BTI® Implants Internal BT	Hex ø 1,22	4.1				
			15.079				non-engaging
			16.079				engaging

AURUMBase®

Angulated screw channel solution. A 25° entry angle, a 360° rotation and a short shaft makes the AURUMBase® the perfect solution for those complicated cases. The Torx® screw combined with the Torx® ball screwdriver allows for maximum torque to be applied at a 25° angle.



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
				36.041	36.042		non-engaging
				37.041	37.042		engaging
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
			36.004	36.005	36.006		non-engaging
			37.004	37.005	37.006		engaging
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
			36.001	36.002			non-engaging
			36.001/2	36.002/2			non-engaging
			36.001/3	36.002/3			non-engaging
			37.001	37.002			engaging
			37.001/2	37.002/2			engaging
			37.001/3	37.002/3			engaging
	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			
			36.007	36.054			non-engaging
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
			36.017	36.018	36.019		
			37.017	37.018	37.019		
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
			36.011	36.012	36.013		
			36.011/2	36.012/2			
			36.011/3	36.012/3			
			37.011	37.012	37.013		
			37.011/2	37.012/2			
			37.011/3	37.012/3			

AURUMBase®


Zimmer® 3i® Biomet®	3i Certain® Implants Internal Hex “Click”	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)	
			36.014	36.015	36.016	non-engaging
			37.014	37.015	37.016	engaging
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.6)	EV (4.2)	EV (4.8)
				36.059	36.060	36.061
				37.059	37.060	37.061
	Astra Tech Osseospeed™ Implants Internal Hex Conic	On Abutment	36.058			
		Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	
		Std. Single		36.059	36.060	non-engaging
		Std. 10 pack		37.059	37.060	engaging
		TIN Single				
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5
		Std. 10 pack	36.038	36.039	36.040	non-engaging
		TIN Single	37.038	37.039	37.040	engaging
		Torx®				
	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
				36.017	36.018	non-engaging
		TIN Single		37.017	37.018	engaging
		Torx®				
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
				36.046	36.047	non-engaging
				37.046	37.047	engaging
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP	
			36.043	36.044		non-engaging
			37.043	37.044		engaging
	Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)	RB/WB	
		GH 1,5mm	36.062	36.063	36.062	non-engaging
		GH 1,5mm	37.062	37.063	37.062	engaging
	Neodent Grand Morse® NEO GM	NEO	GM			
			36.035/3			non-engaging
			37.035/3			engaging

Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0	
				36.065	36.066		non-engaging
				37.065	37.066		engaging
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)		
			36.074	36.075	36.076		non-engaging
			37.074	37.075	37.076		engaging
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
				36.087	36.088	36.89	non-engaging
				37.087	37.088	37.089	engaging
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD			
			36.072	36.073			non-engaging
			37.072	37.073			engaging
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm					
				36.057			non-engaging
				37.057			engaging
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm					
				36.080			non-engaging
				36.180			non-engaging
				37.080			engaging
				37.180			engaging
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28					
				36.057			non-engaging
				37.057			engaging
Phibo®	Phibo® TSH® External Hex PH		NP	SP	WP		
			36.030	36.031	36.032		non-engaging
			36.030/2	36.031/2	36.032/2		non-engaging
			36.030/3	36.031/3	36.032/3		non-engaging
			37.030	37.031	37.032		engaging
			37.030/2	37.031/2	37.032/2		engaging
			37.030/3	37.031/3	37.032/3		engaging

ELLIPTIBase®

Narrow interdental spaces don't represent a challenge anymore. The ELLIPTIBase®, Has been specially developed for narrow implants (3.0-3.5). The elliptical shape achieves almost the same cementation area as a standard Ti-base with the possibility of achieving a more aesthetic result with the 20° entry angle.

- ZrN coating
- SelectGrip® enhanced surface
- Corner/flat versionsm




Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip Corner	3.0	NP (3.5)	RP (4.3)	WP (5.5)
		Flat	56.068 57.068			
Zimmer® 3i Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm Corner	NP (3.5)	RP (4.5)	WP (5.7)	
		Flat	56.017 57.017			
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm Corner	3.0	EV (3.6)	EV (4.2)	EV (4.8)
		Flat	56.406 57.406			
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx® Corner	NNC (3.5)	RN (4.8)	WN (6.5)	
		Flat	56.070 57.070			
	Straumann® Bone Level Conical BL	Torx® Corner	NP (3.3)	RP (4.1)	NP/RP	
		Flat	56.043 57.043			
Henry Schein® Group	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm Corner	3.0	NP (3.5)	RP (4.5)	WP (5.7)
		Flat	56.099 57.099			

C-Base®

The iconic design of the Sirona Ti-base with a twist. Engaging and non-engaging versions, different gingival heights, gold anodised and SelectGrip® surface treatment create the perfect upgrade.

-ExoCad and 3Shape libraries available for download at www.dessdental.com



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
		GH 0,8mm		58.041	58.042		non-engaging
		GH 3,0mm		58.041/3	58.042/3		non-engaging
		GH 0,8mm		59.041	59.042		engaging
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		GH 0,4mm	58.004		58.006	58.097	non-engaging
		GH 0,5mm		58.005			non-engaging
		GH 1,4mm	58.004/2	58.005/2	58.006/2		non-engaging
		GH 3,0mm	58.004/3	58.005/3	58.006/3		non-engaging
		GH 0,4mm	59.004		59.006	59.097	engaging
		GH 0,5mm		59.005			engaging
		GH 1,4mm	59.004/2	59.005/2	59.006/2		engaging
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
		GH 1,5mm	58.001/2	58.002/2			non-engaging
		GH 3,0mm	58.001/3	58.002/3			non-engaging
		GH 0,6mm	59.001	59.002			engaging
		GH 1,5mm	59.001/2	59.002/2			engaging
		GH 3,0mm	59.001/3	59.002/3			engaging
Zimmer® 3i Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
		GH 0,5mm	58.017	58.018	58.019		non-engaging
		GH 1,5mm	58.017/2	58.018/2	58.019/2		non-engaging
		GH 3,0mm	58.017/3	58.018/3	58.019/3		non-engaging
		GH 0,5mm	59.017	59.018	59.019		engaging
		GH 1,5mm	59.017/2	59.018/2	59.019/2		engaging
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
		GH 1,5mm	58.011/2	58.012/2	58.013/2		non-engaging
		GH 3,0mm	58.011/3	58.012/3	58.013/3		non-engaging
		GH 0,6mm	59.011	59.012	59.013		engaging
		GH 1,5mm	59.011/2	59.012/2	59.013/2		engaging
		GH 3,0mm	59.011/3	59.012/3	59.013/3		engaging
	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		
		GH 0,4mm	58.014	58.015			non-engaging
		GH 0,5mm			58.016		non-engaging
		GH 1,5mm	58.014/2	58.015/2	58.016/2		non-engaging
		GH 3,0mm	58.014/3	58.015/3	58.016/3		non-engaging
		GH 0,4mm	59.014	59.015			engaging
		GH 0,5mm			59.016		engaging
		GH 1,5mm	59.014/2	59.015/2	59.016/2		engaging
		GH 3,0mm	59.014/3	59.015/3	59.016/3		engaging

Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.6)	EV (4.2)	EV (4.8)	
		GH 0,8mm	58.059	58.060	58.061	non-engaging
		GH 3,0mm	58.059/3	58.060/3	58.061/3	non-engaging
		GH 0,8mm	59.059	59.060	59.061	engaging
		GH 3,0mm	59.059/3	59.060/3	59.061/3	engaging
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	
		GH 0,8mm		58.024		non-engaging
		GH 0,9mm			58.025	non-engaging
		GH 3,0mm		58.024/3	58.025/3	non-engaging
		GH 0,8mm		59.024		engaging
		GH 0,9mm			59.025	engaging
		GH 3,0mm		59.024/3	59.025/3	engaging
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5
		GH 0,3mm	58.038	58.039		non-engaging
		GH 0,6mm			58.040	58.071
		GH 1,5mm	58.038/2	58.039/2	58.040/2	58.071/2
		GH 3,0mm	58.038/3	58.039/3	58.040/3	58.071/3
		GH 0,3mm	59.038	59.039		engaging
		GH 0,6mm			59.040	59.071
		GH 1,5mm	59.038/2	59.039/2	59.040/2	59.071/2
		GH 3,0mm	59.038/3	59.039/3	59.040/3	59.071/3
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm				
		GH 1,4mm			58.050	non-engaging
		GH 2,5mm			58.050/2	non-engaging
		GH 1,4mm			59.050	engaging
		GH 2,5mm			59.050/2	engaging
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
		GH 1,0mm		58.009		non-engaging
		GH 1,2mm			58.045	non-engaging
		GH 1,0mm		59.009		engaging
		GH 1,2mm			59.045	engaging
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP	
		GH 0,7mm	58.043	58.044		non-engaging
		GH 3,0mm	58.043/3	58.044/3		non-engaging
		GH 0,7mm	59.043	59.044		engaging
		GH 3,0mm	59.043/3	59.044/3		engaging
Henry Schein® Group	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)	
		GH 1,0mm	58.074	58.075	58.076	non-engaging
		GH 2,1mm	58.074/2	58.075/2	58.076/2	non-engaging
		GH 1,0mm	59.074	59.075	59.076	engaging
		GH 2,1mm	59.074/2	59.075/2	59.076/2	engaging
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
		GH 0,5mm		58.087	58.088	58.089
		GH 3,0mm		58.087/3	58.088/3	58.089/3
		GH 0,5mm		59.087	59.088	59.089
		GH 3,0mm		59.087/3	59.088/3	59.089/3
						engaging
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD		
		GH 1,0mm	58.072	58.073		non-engaging
		GH 3,0mm	58.072/3	58.073/3		non-engaging
		GH 1,0mm	59.072	59.073		engaging
		GH 3,0mm	59.072/3	59.073/3		engaging

Abutments

Available in 0°, 15° and 25° with SelectGrip® enhanced surface treatment, the DESS abutments offer an excellent alternative.

- Include final screw 3i



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)
		0°		13.041	13.042	
		15°		22.041	22.042	
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0
		0°	13.004	13.005	13.006	
		15°	22.004	22.005	22.006	
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)	
		0°	13.001	13.002	13.003	
		15°	22.001	22.002	22.003	
	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)	
		0°	13.017	13.018	13.019	
		15°	22.017	22.018	22.019	
Zimmer® 3i® Biomet®	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)	
		0°	13.011	13.012	13.013	
		15°	22.011	22.012	22.013	
	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)	
		0°	13.014	13.015	13.016	
		15°	22.014	22.015	22.016	
		25°		23.015		
Dentsply® Sirona®	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5
		0°	13.038	13.039	13.040	
		15°	22.038	22.039	22.040	
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm				
		GH 0,7mm 0°				13.047
		GH 2,5mm 0°				13.250
		GH 1,5mm 15°				22.048
		GH 2,5mm 15°				22.250
		GH 1,5mm 25°				23.048
		GH 2,5mm 25°				23.250
	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
		0°		13.095	13.096	
		15°		22.095	22.096	
		25°		23.095	23.096	

Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
		0°		13.009*		*inc.captive screw
		15°		22.009*		*inc.captive screw
	Straumann® Bone Level Conical BL	25°		23.009*		*inc.captive screw
		Torx®	NP (3.3)	RP (4.1)	NP/RP	
		0°	13.043	13.044		
Henry Schein® Group	BioHorizons® Implants External Hex BH	15°	22.043	22.044		
		Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)	
		0°	13.020	13.021	13.022	
			22.020	22.021		

Screws

Our Screws have the same geometry as the one’s of the implant manufacturer. Having the same footprint for the screwdriver tip, seat angle, diameter, length, metric, above others, make them identical to the one’s supplied by the implant manufacturer.

- Manufactured in Titanium grade 5 ELI.
- Available in Standard, DLC*, TIN.
- Single unit or 10 packs.



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
		Std. Single	19.068	19.027	19.030	19.030	
		Std. 10 pack		19.027-P10	19.030-P10	19.030-P10	
		DLC	19.268	19.227	19.230	19.230	
		TIN Single	19.668	19.627	19.630	19.630	
		Torx®	19.468	19.441	19.442		AURUMBase® & ELLIPTIBase®
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		Std. Single	19.004	19.005	19.005	19.005	
		10 pack	19.004-P10	19.005-P10	19.005-P10	19.005-P10	
		DLC	19.204	19.205	19.205	19.205	
		TIN Single	19.604	19.605	19.605	19.605	
		Torx®	19.404	19.405	19.405		AURUMBase® & ELLIPTIBase®
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
		Std. Single	19.027	19.028	19.029		
		10 pack	19.027-P10	19.028-P10	19.029-P10		
		DLC	19.227	19.228	19.229		
		TIN Single	19.627	19.628	19.629		
		Torx®	19.401	19.402			AURUMBase® & ELLIPTIBase®
	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			
		Std. Single	19.006	19.031			
		Std. 10 pack	19.006-P10	19.031-P10			
		DLC	19.206	19.231			
		TIN Single	19.606	19.631			
		Torx®	19.306	19.331			AURUMBase® & ELLIPTIBase®
		Hex 1,20mm	19.055				
		Hex 1,22mm	19.065				
		Hex 1,27mm	19.110				

Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
		Std. Single	19.068	19.027	19.030	19.030	
		Std. 10 pack		19.027-P10	19.030-P10	19.030-P10	
		DLC	19.268	19.227	19.230	19.230	
		TIN Single	19.668	19.627	19.630	19.630	
		Torx®	19.468	19.441	19.442		AURUMBase® & ELLIPTIBase®
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		Std. Single	19.004	19.005	19.005	19.005	
		10 pack	19.004-P10	19.005-P10	19.005-P10	19.005-P10	
		DLC	19.204	19.205	19.205	19.205	
		TIN Single	19.604	19.605	19.605	19.605	
		Torx®	19.404	19.405	19.405		AURUMBase® & ELLIPTIBase®
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
		Std. Single	19.027	19.028	19.029		
		10 pack	19.027-P10	19.028-P10	19.029-P10		
		DLC	19.227	19.228	19.229		
		TIN Single	19.627	19.628	19.629		
		Torx®	19.401	19.402			AURUMBase® & ELLIPTIBase®
	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			
		Std. Single	19.006	19.031			
		Std. 10 pack	19.006-P10	19.031-P10			
		DLC	19.206	19.231			
		TIN Single	19.606	19.631			
		Torx®	19.306	19.331			AURUMBase® & ELLIPTIBase®
		Hex 1,20mm	19.055				
		Hex 1,22mm	19.065				
		Hex 1,27mm	19.110				

Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)	
		Std. Single	19.012	19.012	19.012	
		10 pack	19.012-P10	19.012-P10	19.012-P10	
		Torx®	19.417	19.417	19.417	AURUMBase® & ELLIPTIBase®
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)	
		Std. Single	19.034	19.034	19.034	
		10 pack	19.034-P10	19.034-P10	19.034-P10	
		TIN Single	19.634	19.634	19.634	
		Torx®	19.402	19.402	19.402	AURUMBase® & ELLIPTIBase®
	3i Certain® Implants Internal Hex “Click”	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)	
		Std. Single	19.011	19.011	19.011	
		10 pack	19.011-P10	19.011-P10	19.011-P10	
		TIN Single	19.611	19.611	19.611	
		Torx®	19.438	19.438	19.438	AURUMBase® & ELLIPTIBase®

Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.6)	EV (4.2)	EV (4.8)	
		Std. Single	19.059	19.060	19.061	
		TIN Single	19.659	19.659	19.659	
		Torx®	19.441	19.460	19.461	AURUMBase® & ELLIPTIBase®
		On Uniabutment®				
		Std. Single			19.058	
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	
		Std. Single	19.015	19.016	19.017	
		Std. 10 pack		19.016-P10	19.017-P10	
		TIN Single	19.615	19.616	19.617	
		Torx®		19.441	19.425	AURUMBase® & ELLIPTIBase®
		On Uniabutment®				
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5
		Std. Single	19.024	19.024	19.024	19.024
		Std. 10 pack	19.024-P10	19.024-P10	19.024-P10	19.024-P10
		TIN Single	19.624	19.624	19.624	19.624
		Torx®	19.438	19.438	19.438	AURUMBase® & ELLIPTIBase®
		Hex ø 1,00mm				
	Dentsply Ankylos® C/X Implants Internal ANK	Std. Single 1,5mm		19.048		
		Std. 10 pack 1,5mm		19.048-P10		
		TIN Single 1,5mm		19.648		
		Std. Single 2,5mm		19.050		
		C-Base Std Screw		19.702		
		Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
	MIS® Seven Implants Internal Hex MI	Std. Single	19.056	19.081	19.081	19.110
						19.085
		TIN Single	19.656	19.681	19.681	19.685
		Torx®		19.417	19.417	
	MIS® C1 Implants Conic MC	Hex ø 1,27mm	NP (yellow)	SP (purple)	WP (green)	
		Std. Single	19.081	19.056	19.081	
		TIN Single	19.681	19.656	19.681	
	MIS® V3 Implants Conic MV	Hex ø 1,27mm	NP (blue)	SP (purple)		
		Std. Single	19.056	19.081		
		TIN Single	19.656	19.681		
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
		Std. Single	19.070	19.026	19.026	
		10 pack		19.026-P10	19.026-P10	
		TIN Single	19.670	19.626	19.626	
		Torx®	19.441	19.446	19.446	
		On SynOcta®				
		Single		19.009	19.009	
		10 Pack		19.009-P10	19.009-P10	
		TIN Single		19.609	19.609	

Straumann® Group	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP	
		Std. Single	19.032	19.033	19.037	
		10 pack	19.032-P10	19.033-P10	19.037-P10	
		TIN Single	19.632	19.633		
		Torx®	19.443	19.443		
		On SRA®				
	Straumann® BLX Conical BLX	Std. Single			19.098	
		10 pack			19.098-P10	
		TIN Single			19.698	
		Torx®	RB (3.5)	WB (4.5)	RB/WB	
		Std. Single	19.054	19.054	19.054	
		10 pack	19.700	19.700	19.700	
Henry Schein® Group	Neodent Grand Morse® NEO GM	NEO	GM			
		Std. Single	19.051			
	Anthogyr Axiom® Implants Conic Axio	Hex ø 1,18mm				
		Std. Single			19.077	
		TIN Single			19.677	
		On Abutment				
	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0
		Std. Single	19.044	19.044	19.044	19.045
		TIN Single	19.644	19.644	19.644	19.645
		Torx®		19.444	19.444	
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)	
		Std. Single	19.074	19.074	19.076	
		TIN Single	19.674	19.674	19.676	
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
		Std. Single	19.043	19.043	19.043	19.043
		TIN Single	19.643	19.643	19.643	19.643
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD		
		Std. Single	19.072	19.073		
		DLC	19.272	19.273		
		TIN	19.672	19.673		
		Torx®	19.444	19.461		AURUMBase® & ELLIPTIBase®
	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm				
		Std. Single	19.057			
		TIN	19.657			
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm				
		Std. Single	19.080			
		TIN Single	19.680			
		Torx® 0,9mm	19.444			AURUMBase®
		Torx® 2,2mm	19.443			AURUMBase®

GlobalID®	GlobalID® Implants Conic Ink	Hex ø 1,22mm							
		Std. Single	19.086						
		TIN Single	19.686						
		On Abutment							
		Std. Single	19.035						
		TIN Single	19.635						
		Single	19.036				MUA	✓	
		TIN Single	19.636				MUA	✓	
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28							
		Std. Single	19.090						
		TIN Single	19.690						
		Torx®	19.405					AURUMBase®	
		On Abutment							
		Single	19.090						
		TIN Single	19.690						
Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27	3.25/3.75	4.1	4.5	5.5			
		Std. Single	19.091	19.091	19.091	19.091			
		TIN Single	19.691	19.691	19.691	19.691			
		Torx®	19.405					AURUMBase®	
		On Abutment							
		Single	19.090						
		TIN Single	19.690						
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22							
		Std. Single	19.078						
		TIN Single	19.678						
		On Abutment		NP	RP				
		Std. Single		19.063	19.064				
BTI®	BTI® Implants Internal BT	Hex ø 1,22	4.1						
		Std. Single	19.079						
		TIN Single	19.679						
		On Abutment							
		Std. Single		19.065					
		TIN Single		19.665					

Pre-Milled Blanks

DESS pre-milled blanks are available in titanium Grade 5 ELI (ø 10mm - ø 14mm) and Micro Melt (CrCo ø 14mm). With 28 different compatibilities, it is the most comprehensive range in the market.

- Total millable length of 20 mm
- Gingival transition only 0,3 mm to maximize the design freedom
- Includes final screw



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip		3.0	NP (3.5)	RP (4.3)	WP (5.5)	
		ø 10mm	Ti	61.068	61.041	61.042	61.069	
		ø 14mm	Ti		62.041	62.042	62.069	
		ø 14mm	CoCr	71.068	71.041	71.042	71.069	
	Nobel Replace Select™ Trilobe	UniGrip		NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		ø 10mm	Ti	61.004	61.005	61.006	61.097	
		ø 14mm	Ti		62.005	62.006	62.097	
		ø 14mm	CoCr	71.004	71.005	71.006	71.097	
	Nobel Brånemark® External Hex Universal	UniGrip		NP (3.5)	RP (4.1)	WP (5.0)		
		ø 10mm	Ti	61.001	61.002	61.003		
ø 14mm		CoCr	71.001	71.002	71.003			
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm		NP (3.5)	RP (4.5)	WP (5.7)		
		ø 10mm	Ti	61.017	61.018	61.019		
		ø 14mm	Ti	62.017	62.018	62.019		
		ø 14mm	CoCr	71.017	71.018	71.019		
	3i Osseotite® External Hex USA	Hex ø 1,27mm		NP (3.4)	RP (4.1)	WP (5.0)		
		ø 10mm	Ti	61.011	61.012	61.013		
		ø 14mm	CoCr	71.011	71.012	71.013		
	3i Certain® Implants Internal Hex “Click”	Hex ø 1,20mm		NP (3.4)	RP (4.1)	WP (5.0)		
		ø 10mm	Ti	61.014	61.015	61.016		
		ø 14mm	Ti		62.015	62.016		
ø 14mm		CoCr	71.014	71.015	71.016			
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.06)	EV (3.6)	EV (4.2)	EV (4.8)	EV (5.4)
		ø 10mm	Ti	61.406	61.059	61.060	61.060	
		ø 14mm	Ti		62.059	62.060	62.060	62.407
		ø 14mm	CoCr		71.059	71.060	71.060	71.407
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm		NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)		
		ø 10mm	Ti	61.023	61.024	61.025		
		ø 14mm	Ti		62.024	62.035		
		ø 14mm	CoCr		71.024	71.025		
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm		NP (3.4)	RP (3.8)	WP (4.5)	5.5	
		ø 10mm	Ti	61.038	61.039	61.040	61.071	
		ø 14mm	Ti		62.039	62.040	62.071	
		ø 14mm	CoCr	71.038	71.039	71.040	71.071	
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm						
		ø 10mm	Ti	61.038				
ø 14mm		CoCr	71.038					

Pre-Milled Blanks

Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm		NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
		ø 10mm	Ti	61.056	61.095	61.096	
		ø 14mm	Ti				
	MIS® C1 Implants Conic MC	Hex ø 1,27mm		NP (yellow)	SP (purple)	WP (green)	
		ø 10mm	Ti	61.081	61.082	61.083	
	MIS® V3 Implants Conic MV	Hex ø 1,27mm		NP (blue)	SP (purple)		
		ø 10mm	Ti	61.084	61.082		
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®		NNC (3.5)	RN (4.8)	WN (6.5)	
		ø 10mm	Ti	61.070	61.009	61.045	
		ø 14mm	Ti		62.009	62.045	
	Straumann® Bone Level Conical BL	Torx®		NP (3.3)	RP (4.1)	NP/RP	
		ø 10mm	Ti	61.043	61.044		
		ø 14mm	Ti	62.043	62.044		
	Straumann® BLX Conical BLX	Torx®		RB (3.5)	WB (4.5)	RB/WB	
		ø 10mm	Ti	61.062	61.062	61.062	
		ø 14mm	Ti		62.063		
	Neodent Grand Morse® NEO GM	NEO		GM			
		ø 10mm	Ti	61.035			
		ø 14mm	Ti	62.035			
	Anthogyr Axiom® Implants Conic Axio	Hex ø 1,18mm					
		ø 10mm	Ti	61.077			
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm		3.3	3.8	4.3	5.0
		ø 10mm	Ti	61.064	61.065	61.066	61.067
		ø 14mm	Ti		62.065	62.066	62.067
	Conelog® Implants Internal Con	Hex ø 1,27mm		NP (3.3)	RP (3.8)	WP (5.0)	
		ø 10mm	Ti	61.074	61.075	61.076	
		ø 14mm	Ti	62.074	62.075	62.076	
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm		3.0	NP (3.5)	RP (4.5)	WP (5.7)
		ø 10mm	Ti	61.099	61.087	61.088	61.089
		ø 14mm	Ti		62.087	62.088	62.089
	BioHorizons® Implants External External Hex BH	Hex ø 1,27mm		NP (3.7)	RP (4.2)	WP (5.2)	
		ø 10mm	Ti	61.020	61.021	61.022	
		ø 14mm	Ti	62.020	62.021	62.022	

Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm		MINI	STD		
		ø 10mm	Ti	61.072	61.073		
		ø 14mm	Ti	62.072	62.073		
		ø 14mm	CoCr	71.072	71.073		
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm					
		ø 10mm	Ti	61.057			
		ø 14mm	CoCr	71.057			
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm					
		ø 10mm	Ti	61.080			
		ø 14mm	Ti	62.080			
		ø 14mm	CoCr	71.080			
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28					
		ø 10mm	Ti	61.080			
		ø 14mm	Ti	62.080			
		ø 14mm	CoCr	71.080			
Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27		3.25/3.75	4.1	4.5	5.5
		ø 10mm	Ti	61.091	61.092	61.093	61.094
		ø 14mm	Ti	62.094			
		ø 14mm	CoCr	71.094			
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22					
		ø 10mm	Ti	61.078			
		ø 14mm	Ti				
		ø 14mm	CoCr	71.078			
Blocks	Blocks						
		Ti		63.020		□15mm	L 20mm
		Ti		63.040		□15mm	L 40mm
		CoCr		73.020		□15mm	L 20mm
		CoCr		73.040		□15mm	L 40mm

CoCr bases

Manufactured in micromelt CoCr alloy with a sanitary POM sleeve, the DESS CoCr bases are perfect for casting techniques. With a higher melting point (1369°C), they will maintain precision of shape in the fusion with other metals.

- Can be used in bone-level restorations or from abutment level (Multiunit)
- Engaging and non-engaging versions available
- Final screw includedmet



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
				30.041	30.042		non-engaging
				31.041	31.042		engaging
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
			30.004	30.005	30.006		non-engaging
			31.004	31.005	31.006		engaging
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
			30.001	30.002	30.003		non-engaging
			31.001	31.002	31.003		engaging
				30.102			non-engaging For welding techniques
	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			
			30.007	30.054			non-engaging
		Hex 1,20mm	30.007H120				non-engaging
		Hex 1,22mm	30.007H122				non-engaging
		Hex 1,27mm	30.007H127				non-engaging
		Torx® Ball	30.007TB				non-engaging
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
			30.017	30.018	30.019		non-engaging
			31.017	31.018	31.019		engaging
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
			30.011	30.012	30.013		non-engaging
			31.011	31.012	31.013		engaging
	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		
			30.014	30.015	30.016		non-engaging
			31.014	31.015	31.016		engaging

Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.6)	EV (4.2)	EV (4.8)	
				30.059	30.060	30.061	non-engaging
				31.059	31.060	31.061	engaging
	Astra Tech Osseospeed™ Implants Internal Hex Conic	On Uniabutment® Hex ø 1,27mm	30.058	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	
					30.024	30.025	non-engaging
					31.024	31.025	engaging
		On Uniabutment® On Uniabutment®	30.036				non-engaging
			30.037				non-engaging
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5	
			30.038	30.039	30.040		non-engaging
			31.038	31.039	31.040		engaging
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm					
		GH 1,5mm				30.050	non-engaging
		GH 2,5mm				30.250	non-engaging
		GH 1,5mm				31.050	engaging
		GH 2,5mm				31.250	engaging
	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8	
				30.056	30.095	30.096	non-engaging
				31.056	31.095	31.096	engaging
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)		
				30.009	30.045		
				31.009	31.045		non-engaging For welding techniques
				30.109	30.145		non-engaging For welding techniques
		On Abutment		31.109	31.145		engaging
		On Abutment		30.010	30.048		engaging
	Straumann® Bone Level Conical BL			31.010	31.048		non-engaging For welding techniques
				30.110			
		Torx®	NP (3.3)	RP (4.1)	NP/RP		non-engaging
			30.043	30.044			engaging
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0	
				30.065	30.066		non-engaging
				31.065	31.066		engaging
	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD			
			30.072	30.073			non-engaging
			31.072	31.073			engaging
	Dentium® Superline™ & Implantium® Conic Dent	Hex ø 1,28					
					30.090		non-engaging
					31.090		engaging

Castables

Unlike other products on the market, DESS castables are machined, which produces a more accurate connection replica, avoiding imperfections in the casting.

- Engaging and non-engaging
- Single and 10 unit packs
- Single unit includes final screw



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
				28.041	28.042		non-engaging
				11.041-P10	11.042-P10		
				29.041	29.042		engaging
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
			28.004	28.005	28.006		non-engaging
			11.004-P10	11.005-P10			
			29.004	29.005	29.006		engaging
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
			28.001	28.002	28.003		non-engaging
				11.002-P10			
			29.001	29.002	29.003		engaging
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
			28.001	28.002	28.003		non-engaging
				11.002-P10			
			29.001	29.002	29.003		engaging
	3i Osseotite® External Hex USA	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
			28.001	28.002	28.003		non-engaging
				11.002-P10			
			29.001	29.002	29.003		engaging
	3i Certain® Implants Internal Hex "Click"	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
			28.001	28.002	28.003		non-engaging
				11.002-P10			
			29.001	29.002	29.003		engaging

Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.6)	EV (4.2)	EV (4.8)	
		On Uniabutment®	28.059	28.060	28.061	non-engaging
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	
				28.024	28.025	non-engaging
				11.024-P10	11.025-P10	
				29.024	29.025	engaging
				12.024-P10	12.025-P10	
		On Uniabutment®	28.036			non-engaging 20°
	Dentsply Friadent® Xive® Implants Internal Hex FD	On Uniabutment®	11.036-P10			20°
		On Uniabutment®	28.037			non-engaging 45°
		On Uniabutment®	11.037-P10			45°
		Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5
			28.038	28.039	28.040	non-engaging
				11.039-P10		
Straumann® Group	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm				
		GH 1,5mm		28.050		
	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
				28.056	28.095	28.096
					11.095-P10	11.096-P10
				29.056	29.095	29.096
					12.095-P10	12.096-P10
	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
				28.009	28.045	non-engaging
				11.009-P10		
				29.009	29.045	engaging
				12.009-P10		
		On Abutment		28.010	28.046	non-engaging
Henry Schein® Group	Straumann® Bone Level Conical BL	On Abutment		11.010-P10		
		On Abutment		29.010	29.046	engaging
		On Abutment		12.010-P10		
		Torx®	NP (3.3)	RP (4.1)	NP/RP	
			28.043	28.044		non-engaging
			11.043-P10	11.044-P10		
	Camlog® Implants Internal Cam			29.043	29.044	engaging
				12.043-P10	12.044-P10	
		Hex ø 1,27mm	3.3	3.8	4.3	5.0
				28.065	28.066	non-engaging
Dentium®	BioHorizons® Implants External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)	
		Std. Single	28.020	28.021	28.022	non-engaging
		Std. 10 pack		11.021-P10		
		TIN Single	29.020	29.021	29.022	engaging
				12.021-P10		
	Superline™ & Implantium® Conic Dent	Hex ø 1,28				
				28.090		non-engaging
				29.090		engaging

Ti-Base Castables

Manufactured in sanitary POM, the Ti-base sleeves are perfect for casting processes.

- Available in packs of 5
- Engaging and non-engaging versions



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
			32.001-P5	32.001-P5	32.001-P5	32.001-P5	non-engaging
			33.001-P5	33.001-P5	33.001-P5	33.001-P5	engaging
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
			32.001-P5	32.001-P5	32.001-P5	32.001-P5	non-engaging
			33.001-P5	33.001-P5	33.001-P5	33.001-P5	engaging
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
			32.001-P5	32.001-P5	32.003-P5		non-engaging
			33.001-P5	33.001-P5	33.003-P5		engaging
	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			
			32.007-P5	32.054-P5			non-engaging
			32.007/8-P5				engaging
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
			32.001-P5	32.001-P5	32.001-P5		non-engaging
			33.001-P5	33.001-P5	33.001-P5		engaging
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
			32.001-P5	32.001-P5	32.001-P5		non-engaging
			33.001-P5	33.001-P5	33.001-P5		engaging
	3i Certain® Implants Internal Hex “Click”	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		
			32.001-P5	32.001-P5	32.001-P5		non-engaging
			33.001-P5	33.001-P5	33.001-P5		engaging
	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.6)	EV (4.2)	EV (4.8)	
				32.001-P5	32.001-P5	32.001-P5	non-engaging
				33.001-P5	33.001-P5	33.001-P5	engaging
		On Uniabutment®	32.001-P5				non-engaging
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm GH 1,1mm		NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	
				32.001-P5	32.001-P5	32.001-P5	non-engaging
				33.001-P5	33.001-P5	33.001-P5	engaging
		On Uniabutment®	32.001-P5				non-engaging
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5	
			32.001-P5	32.001-P5	32.001-P5	32.001-P5	non-engaging
			33.001-P5	33.001-P5	33.001-P5	33.001-P5	engaging
Dentsply Sirona®	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm				32.001-P5	non-engaging
						33.001-P5	engaging


Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8	
			32.001-P5	32.001-P5	32.001-P5	32.007-P5	non-engaging
			33.001-P5	33.001-P5	33.001-P5		engaging
	MIS® C1 Implants Conic MC	Hex ø 1,27mm	NP (yellow)	SP (purple)	WP (green)	Multi-unit 4.8	
			32.001-P5	32.001-P5	32.001-P5	32.007-P5	non-engaging
			33.001-P5	33.001-P5	33.001-P5		engaging
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)		
			32.001-P5	32.046-P5	32.045-P5		non-engaging
			33.001-P5	33.009-P5	33.045-P5		non-engaging
		On Abutment On Abutment		32.010-P5	32.045-P5		non-engaging
		Torx®		33.010-P5	33.045-P5		engaging
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP		
			32.001-P5	32.001-P5			non-engaging
		On SRA On SRA					engaging
					32.098-P5		non-engaging
	Straumann® BLX Conical BLX	Torx®			33.098-P5		engaging
		Torx®	RB (3.5)	WB (4.5)	RB/WB		
			32.001-P5	32.001-P5	32.001-P5		non-engaging
	Neodent Grand Morse® NEO GM	NEO					engaging
					32.001-P5		non-engaging
		Hex ø 1,18mm			33.001-P5		engaging
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0	
			32.001-P5	32.001-P5	32.001-P5	32.001-P5	non-engaging
			33.001-P5	33.001-P5	33.001-P5	33.001-P5	engaging
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)		
			32.001-P5	32.001-P5	32.001-P5	32.001-P5	non-engaging
			33.001-P5	33.001-P5	33.001-P5	33.001-P5	engaging
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
				32.001-P5	32.001-P5	32.001-P5	non-engaging
				33.001-P5	33.001-P5	33.001-P5	engaging
	BioHorizons® Implants External External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)		
			32.001-P5	32.001-P5	32.001-P5		non-engaging
			33.001-P5	33.001-P5	33.001-P5		engaging

Ti-Base Castables

Hiossen®/Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI		STD		non-engaging engaging
			32.001-P5	32.001-P5			
			33.001-P5	33.001-P5			
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm			32.001-P5		non-engaging
					33.001-P5		engaging
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm			32.001-P5		non-engaging
					33.001-P5		engaging
GlobalD®	GlobalD® Implants Conic Ink				32.001-P5		non-engaging
					33.001-P5		engaging
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28			32.001-P5		non-engaging
					33.001-P5		engaging
Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27	3.25/3.75	4.1	4.5	5.5	non-engaging
			32.001-P5	32.001-P5	32.001-P5	32.001-P5	engaging
			33.001-P5	33.001-P5	33.001-P5	33.001-P5	
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22		NP	RP		non-engaging
			32.001-P5				engaging
			33.001-P5				
BTI®	BTI® Implants Internal BT	Hex ø 1,22	4.1				non-engaging
			32.001-P5				engaging
			33.001-P5				

AURUMBase® Castables

Machine manufactured in sanitary POM and available in 0°, 10°, 20°, the AURUMBase® castables, will facilitate the creation of a precise angulated channel when used on traditional workflows.
- Available in packs of 5 units



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)
		0°		33.100-P5	33.100-P5	
		10°		33.101-P5	33.101-P5	
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0
		0°	33.100-P5	33.100-P5	33.100-P5	
		10°	33.101-P5	33.101-P5	33.101-P5	
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)	
		0°	33.100-P5	33.100-P5		
		10°	33.101-P5	33.101-P5		
	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)		
		0°	33.700-P5	33.700-P5		
		10°	33.701-P5	33.701-P5		
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)	
		0°	33.700-P5	33.700-P5		
		10°	33.701-P5	33.701-P5		
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)	
		0°	33.100-P5	33.100-P5	33.100-P5	
		10°	33.101-P5	33.101-P5	33.101-P5	
	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)	
		0°	33.100-P5	33.100-P5	33.100-P5	
		10°	33.101-P5	33.101-P5	33.101-P5	
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.6)	EV (4.2)	EV (4.8)	
		0°	33.100-P5	33.100-P5	33.100-P5	
		10°	33.101-P5	33.101-P5	33.101-P5	
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	
		0°		33.100-P5	33.100-P5	
		10°		33.101-P5	33.101-P5	
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5
		0°	33.100-P5	33.100-P5		
		10°	33.101-P5	33.101-P5		
	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
		0°		33.100-P5	33.100-P5	
		10°		33.101-P5	33.101-P5	

AURUMBase® Castables

Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
		0°		33.200-P5	33.300-P5	
		10°		33.201-P5	33.301-P5	
		20°		33.202-P5	33.302-P5	
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP	
		0°	33.100-P5	33.100-P5		
		10°	33.101-P5	33.101-P5		
		20°	33.102-P5	33.102-P5		
	Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)	RB/WB	
		0°	33.100-P5	33.100-P5	33.100-P5	
		10°	33.101-P5	33.101-P5	33.101-P5	
		20°	33.102-P5	33.102-P5	33.102-P5	
Neodent Grand Morse® NEO GM	NEO	GM				
	0°	33.100-P5				
	10°	33.101-P5				
	20°	33.102-P5				
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0
		0°		33.100-P5	33.100-P5	
		10°		33.101-P5	33.101-P5	
		20°		33.102-P5	33.102-P5	
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)	
		0°	33.100-P5	33.100-P5	33.100-P5	
		10°	33.101-P5	33.101-P5	33.101-P5	
		20°	33.102-P5	33.102-P5	33.102-P5	
	BioHorizons® Implants Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
		0°		33.100-P5	33.100-P5	
		10°		33.101-P5	33.101-P5	
		20°		33.102-P5	33.102-P5	
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD		
		0°	33.100-P5	33.100-P5		
		10°	33.101-P5	33.101-P5		
		20°	33.102-P5	33.102-P5		
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm				
		0°	33.100-P5			
		10°	33.101-P5			
		20°	33.102-P5			
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm				
		0°	33.100-P5			
		10°	33.101-P5			
		20°	33.102-P5			
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28				
		0°	33.100-P5			
		10°	33.101-P5			
		20°	33.102-P5			
Phibo®	Phibo® TSH® External Hex PH		NP	SP	WP	
		0°	33.100-P5			
		10°	33.101-P5			
		20°	33.102-P5			





Holders

Holders		For standard 98mm disc machines	Any milling machine with a 98 mm disc holder
	HLDSTD		
		For Imes-Icore® machine	Coritec 245i touch / Coritec 250i touch Coritec 350i touch / Coritec 350i touch Loader Coritec 450i / Coritec 550i Coritec 650i / Coritec 650i Loader
	HLDIMES		
		For VHF®	VHF® S1 VHF® S2 VHF® J5
	HLDVHF		
		For VHF®	VHF® R5
	HLDVHFR5		
		For YENADENT®	D6 D15 D43N
	HLDYENA		
		For VHF®	VHF® N4
	HLDVHFN4		
		For VHF®	VHF® Z4
	HLDVHFZ4		
		For machines with support / changer 3R	D14P DC40 D43N
	HLD3R		

Tools


Our screws have the same geometry as the implant manufacturer's ones. Having the same footprint for the screwdriver tip, seat angle, diameter, length, metric, and more makes them identical to the ones supplied by the implant manufacturer.

- Manufactured in titanium Grade 5 ELI
- Available in Standard, DLC*, TIN
- Single unit or packs of 10

Screwdrivers	HEXAGON	18mm	20mm	24mm	30mm	35mm		
		0,90 mm	DTC20H090	DTC24H090				
		1,00 mm	DTC20H100	DTC24H100		DTC35H100		
		1,20 mm	DTC18H120	DTC20H120	DTC24H120	DTC30H120		DTC35H120
		1,22 mm	DTC20H122	DTC24H122		DTC35H122		
		1,25 mm	DTC20H125	DTC24H125				
		1,27 mm	DTC20H127	DTC24H127	DTC30H127	DTC35H127		
		1,30 mm	DTC20H130	DTC24H130				
		1,50 mm	DTC20H150	DTC24H130				
		1,70 mm	DTC20H170	DTC24H170				
1,75 mm	DTC20H175	DTC24H175						
SQ		DT20S130	DT24S130					
Torx® 06		DT20T06	DT24T06	DT30T06	DT35T06			
		DT20TC06	DT24TC06	DT30TC06	DT35TC06			
UG		DT20TB06	DT24TB06	DT30TB06	DT35TB06			
NEO			DT24NG	DT30NG	DT35NG			
Screwdriver handles	Premium laboratory handle		PHLABGOLD					
			PHLABBLACK					
	Dentist handle	7 mm	MDTMAN700			DT35T06		
	9 mm	MDTMAN900			DT35TC06			

Tools

Special tools	Multi-unit® RP	DTMURP	
	Multi-unit® WP	DTMUWP	
	UniAbutment®	DTUNI	
	Octagon	DTSYN	
	Torque wrench adaptor	TWADAP	
	DESSLoc®	DTLOC	
	Adaptor from 6 mm to straumann® torque wrench to square conection	DTSTSQ	
		DTSTWH	
	Screwdriver tips adaptor fro straumn® torque wrench	DTACTLOC	
		DTSTRLOC	
	Starumann® torque wrench	DTSTCA	

Special tools	Nobel Biocare™ adaptor from 8mm to square connection	DTNBWH	
	Nobel Biocare™ adapter from 8mm to square connection	DTNBSQ	
	Screwdriver tips adaptor from Nobel Biocare™ torque wrench	DTNBWH	
	Nobel Biocare™ contra-angles adapter	DTNBCA	
	Handpiece adaptor for square connection	DTANGSQ	
	Manual driver for Intraoral Scan abutment	DTPEI	
	Manual tool for intraoral Scan abutment	DTPEIPEEK-P2	



DESS HQ



6.200m² Barberà del Vallès, Barcelona Spanien



All DESS products are manufactured according to ISO 9001 and ISO 13485 with CE marking for all products Class IIb

All DESS® AURUMBase® screws and screwdrivers are based on the standard ISO 10664

3i® and 3i Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations, USA

ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XiVE®, FRIADENT®, ASTRATECH implant system™ EV and UniAbutment®

Cerec® are registered trademarks of DENTSPLY SIRONA Implants

BioHorizons® is a registered trademark of BioHorizons, INC.

Brånemark System®, NobelReplace® and Replace Select™, NobelActive™, NobelReplace® CC, NobelSpeedy™ and Multi-unit® are registered trademarks of Nobel Biocare AB,

Neodent®, Grand Morse®, Straumann® and synOcta® are registered trademarks of Straumann Holding AG, Switzerland

MIS® is a registered trademark of MIS Implant Technologies Ltd.

Zimmer Screw-Vent® and SwissPlus® are registered trademarks of Zimmer Dental INC, USA

Dyna® is a registered trademark of Dyna Dental Engineering BV

Locator® is a registered trademark of ZEST IP Holdings

Camlog® and Conelog® are registered trademarks of CAMLOG Biotechnologies AG

Megagen® is a registered trademark of MEGAGEN IMPLANT CO. LTD

TORX® is a trademark of ACUMENT INTELLECTUAL PROPERTIES, LLC

Bego Semados® is a registered trademark of BEGO GmbH & Co. KG

Anthogyr® is a registered trademark of Anthogyr SAS

Medentis ICX® is a registered trademark of medentis medical GmbH

Osstem®, Hiossen® is a registered trademark of OSSTEM IMPLANT CO. LTD

Dentium® is a trademark of Dentium Co. Ltd.

GOBALD® is a trademark of Global D SASU.

Biotech Kontakt® is a trademark Biotech Dental SAS.

BTI® is a trademark BTI Biotechnology Institute.

PHIBO® is a trademark PHIBO Dental Solutions S.L.

DESS products can be used only by licensed dentists, physicians and labs.
Please refer ever to the instructions of use available on the website
www.dessdental.com
or supplied with the product.

Some products may not be available in your country. Some of the products are cleared for US and Canada.
Please contact your local distributor for more details.



